


Historic, Archive Document

Do not assume content reflects current scientific knowledge, policies, or practices.

IMMUNOELECTROOSMOPHORESIS - IEOP
AGAR GEL DIFFUSION PRECIPITATION - AGDP
TESTS

aSF977
.A4D5
v.3



U. S. DEPT. OF AGRICULTURE
NATIONAL AGRICULTURAL LIBRARY

JUL 14 1981


CATALOGING = PREP

PLUM ISLAND ANIMAL DISEASE CENTER
POST OFFICE BOX 848
GREENPORT, NEW YORK 11944

AD-33 Bookplate
(1-63)

NATIONAL

**A
G
R
I
C
U
L
T
U
R
A
L**



LIBRARY

Introduction

Control and eradication of African swine fever (ASF) depend upon detection and elimination of exposed animals. The hemadsorption (HAd) test*, as well as immunofluorescence are valuable in the diagnosis of acute ASF. Where the disease has become enzootic in domestic swine, subacute and chronic infections are often encountered. Chronic infections stimulate complement fixing (CF) and AGDP antibodies against ASF virus, but detection by their use requires at least 12 hours. The IEOP test is more rapid and sensitive than either the CF or AGDP tests in detecting antibodies against ASF virus. (Pan, I. C., De Boer, C. J., and W. R. Hess.

*Please see Microfiche PIL-M-1.

African Swine Fever: Application of
Immuno-electro-osmophoresis for Detection
of Antibody. Canadian Journal of
Comparative Medicine, 36:309-316, 1972.)

Since it is seldom possible to
establish a diagnosis of ASF by clinical
means alone, laboratory tests are re-
quired to differentiate the disease from
hog cholera and other diseases. This
microfiche describes in detail how the
IEOP, AGDP and FA tests for ASF are per-
formed. Laboratory personnel with
experience in cell culture (CC) should
be able to carry out the diagnostic
procedures described.

1951
African Swine Fever: Application of
Immunoelectrophoresis for Detection
of Antibody. Canadian Journal of
Comparative Medicine, 20:309-316, 1952.
Since it is not possible to
establish a diagnosis of ASF by clinical
means alone, laboratory tests are re-
quired to differentiate the disease from
hog cholera and other diseases. This
microscopic description is useful for the
identification of the virus. The virus
is found in laboratory personnel with
experience in cell culture (CC) should
be able to carry out the diagnosis
readily and accurately.

THE IEOP TEST

The next frame illustrates graphically the basic steps in the IEOP test. The test is fundamentally simple and is easily performed if the laboratory has good ASF viral antigen and an immunoelectrophoresis apparatus. Microscopic slides are coated with agar in which wells are made. Antigen is placed in wells nearest the cathode, or negative pole and serum in the opposite wells. Direct electric current is applied for about 30 minutes. At the end of that period the test may be read.

THE TEST

The next frame illustrated graphically

only the basic steps in the test. The test is fundamentally simple and is easily performed at the laboratory and good A-B visual antigen and an immune electrophoretic apparatus. Microscopic slides are coated with agar in which wells are made. Antigen is placed in wells against the cathode, or negative pole and serum in the opposite wells. Direct electric current is applied for about 30 minutes. At the end of that period the test may be read.


BASIC STEPS IN THE IMMUNOELECTROOSMOPHORESIS (IEOP) TEST FOR AFRICAN SWINE FEVER (ASF) ANTIBODIES

1. Ordinary microscope slides are coated with agar gel.




Diameter of well= 2.5mm
Interwell distance (center to center) = 13mm


2. Wells are made in the agar.




3. Known ASF antigen (AG) is placed in one of a pair of wells and serum(s) from a suspect animal in the other. (Control antisera are also used.)



4. Direct current is applied to the slide with the negative (-, cathode) pole connected to the side with antigen for 30 minutes.



5. Where the serum is positive, lines appear between the wells.



(Advantages: More than 100 tests can be performed at once. The test is accurate and sensitive requiring small amounts of reagents.)



POSITIVE IEOP TESTS

Wells on the left were filled with antisera ("Lisbon" , "Rome") from swine infected with ASF virus.

Wells on the right were filled with ASF viral antigen. Lines in the agar between the wells indicate positive sera.

TABLE I. The Sensitivity of the IEOP Test by Block Titration of Antigen and Antibody

Ag									
As		u*	1:2	1:4	1:8	1:16	1:32	1:64	1:128
u.....		+	+	+	+	+	-	-	-
1:2.....		+	+	+	+	+	+	-	-
1:4.....		+	+	+	+	+	+	-	-
1:8.....		+	+	+	+	+	+	-	-
1:16.....		+	+	+	+	+	+	+	-
1:32.....		+	+	+	+	+	+	+	-
1:64.....		+	+	+	+	+	+	-	-
1:128.....		+	+	+	+	+	+	-	-
1:256.....		-	+	+	+	+	+	-	-
1:512.....		-	-	+	+	+	+	-	-
1:1024.....		-	-	-	+	+	-	-	-
1:2048.....		-	-	-	-	-	-	-	-

Note: Reagents volume = 5 μ l. Final results recorded after staining

*u = undiluted

TABLE II. The Sensitivity of the CF Test* by Block Titration of Antigen and Antibody.

Ag							
As		1:8	1:16	1:32	1:64	1:128	1:256
1:8.....		+	+	+	+	-	-
1:16.....		+	+	+	+	+	-
1:32.....		+	+	+	+	+	-
1:64.....		+	+	+	+	-	-
1:128.....		+	+	+	-	-	-
1:256.....		+	+	+	-	-	-
1:512.....		+	+	-	-	-	-
1:1024.....		-	-	-	-	-	-
1:2048.....		-	-	-	-	-	-

*The reaction mixtures contained 0.25 ml of serum dilution, 1.75 hemolytic units of guinea pig complement in 0.5 ml and 0.25 ml antigen dilution. Fixation was overnight at 5°C. Sensitized sheep red blood cells (2%) (0.5 ml) were added and the mixtures incubated for one hour at 37°C. Tubes which showed complete inhibition of hemolysis were expressed as +

TABLE III. The Sensitivity of the AGDP Test^a by Block Titration of Antigen and Antibody

Ag				
As	Undiluted	1:2	1:4	1:8
Undiluted	+	+	+	-
1:2.....	+	+	+	-
1:4.....	+	+	+	-
1:8.....	+	+	+	-
1:16.....	+	+	-	-
1:32.....	+	+	-	-
1:64.....	+	+	-	-
1:128.....	+	+	-	-
1:256.....	+	+	-	-
1:512.....	-	-	-	-

^aReagent volume = 30 μ l. Well size was 4 mm in diameter, and the distance of antigen and anti-serum was 4 mm. Results recorded five days incubation at room temperature.

TABLE IV. Detection of Antibody in Pig Sera by Three Different Tests

Infected With	IEOP			CF	AGDP	IEOP (+)^b
	Undiluted^a	Diluted	Staining			AGDP (+) CF (-)
Lisbon	115/120 ^a (95.8%)	119/120 (99.1%)	120/120 (100%)	85/120 (70.8%)	84/120 (70.0%)	5/120 (4.1%)
Tengani	95/95 (100%)	95/95 (100%)	95/95 (100%)	33/67 (47.2%)	90/93 (96.7%)	31/67 (46.2%)
Rome	48/52 (92.3%)	52/52 (100%)	52/52 (100%)	30/52 (57.5%)	48/52 (92.3%)	27/52 (42.3%)
Total	258/267 (96.6%)	266/267 (99.6%)	267/267 (100%)	148/239 (61.9%)	222/265 (83.7%)	63/239 (26.4%)

^aUndiluted antigen was used in screening. Negative sera were retested with diluted antigen

^bSera positive in IEOP and AGDP tests, but negative in the CF test

^cNumerator = number positive; denominator = number tested

Frame A-6 shows the appearance of the IEOP test. Frames A-7 and A-8 summarize typical results of the IEOP test in-comparison with the complement fixation (CF) and AGDP test. An antibody titer of 1:1024 was attained (Table I). This compared favorably to an antibody titer of 1:512 in the CF test (Table II) and 1:256 in the AGDP test (Table III). Using sera from known cases of ASF, the IEOP test revealed a higher number of positives than the others (Table IV).*

*Data in frames are reproduced with the kind permission of the editor of the CANADIAN JOURNAL OF COMPARATIVE MEDICINE, and also Doctors I. C. Pan, C. J. De Boer and W. R. Hess of PIADL.

ANTIGEN REQUIRED FOR IEOP AND AGDP TESTS

It is necessary to have antigen of good quality and high titer for these tests. Laboratories in the USA other than the Plum Island Animal Disease Laboratory (PIADL) are not allowed to possess or make ASF viral antigen. In the USA these tests are performed at PIADL.

Brief but detailed protocols for preparation of the necessary antigen and performing these tests are included for the benefit of laboratories and personnel where the above restrictions do not apply. Laboratory technicians who have had experience in CC and in serologic tests should be able to prepare the antigen and carry out the tests.

ANTHONY ROSSIGNOL FOR THE DEFENSE

It is necessary to have a copy of the
good quality and high quality film for the
best of the film in the FBI office
than the film in the FBI office
of the FBI office and of the FBI office
in the FBI office and of the FBI office
in the FBI office and of the FBI office

END

first and detailed protocol for
the production of the necessary material
performing these tasks and their
the results of these studies and their
where the above mentioned is not
this is a very important task and
to be done and to be done in a
in a very important task and
in a very important task and

Directors of laboratories who wish to have personnel trained in the performance of these tests are invited to write to the Director, USDA, SEA, Plum Island Animal Disease Center, P. O. Box 848, Greenport, New York 11944, USA.

Note: "Mention of a trademark or proprietary product does not constitute a guarantee or warranty of the product by the U.S. Department of Agriculture and does not imply its approval to the exclusion of other products that may also be suitable."

Department of Labor, Bureau of Labor Statistics

It is noted that the Bureau of Labor Statistics has been requested to conduct a study of the effect of the minimum wage law on the cost of living in the District of Columbia. The Bureau has been requested to conduct a study of the effect of the minimum wage law on the cost of living in the District of Columbia.

Washington, D. C. 20540, USA

It is noted that the Bureau of Labor Statistics has been requested to conduct a study of the effect of the minimum wage law on the cost of living in the District of Columbia. The Bureau has been requested to conduct a study of the effect of the minimum wage law on the cost of living in the District of Columbia.

Page 1

Preparation of Antigen

For the IEOP and AGDP Tests

Antigen may be prepared by using a stable line of monkey kidney cells, such as the MS or Vero lines. The MS line of cells was obtained from the Razi Institute of Tehran, Iran and originally came from the National Institute of Animal Diseases in Tokyo, Japan. The Vero line of African green monkey kidney cells was obtained from Dr. A. J. Kniazeff, Naval Biological Laboratory, Berkeley, California in 1967 and has undergone more than 220 passages since then. The monolayer CC are prepared in Povitsky bottles (Corning Glass Works, Corning, N. Y.) (Pan et al., 1972) and the Vero cells in Baxter bottles (Baxter Laboratories, Morton Grove, Illinois). Drums of Baxter bottles on a roller mill are shown in Frame B-6. The recovery

1900

1901

1902

1903

1904

1905

1906

1907

1908

1909

1910

1911

1912

1913

1914

1915

1916

1917

1918

1919

1920

1921

1922

1923

1924

1925

1926

of ASF antigen from the Vero line in Baxter bottles is better than from MS cells in Povitsky bottles. A roller mill is required to make monolayers in Baxter bottles. Other types of roller mills with large bottles may also be used. Both cell lines can also be cultured in Povitsky bottles, plastic CC containers, or in prescription bottles as stationary cultures; amounts of antigen produced by these means may be adequate.

However, since the roller mill method employing Baxter bottles is highly productive and is now employed at the Plum Island Animal Disease Laboratory, it will be described in detail.

1. Preparation of Growth Medium for Vero
Cells

1.1 A sterile 20 liter (L) carboy
(large glass bottle or container) is
filled with 10 L of distilled (D) H₂O.

1.2 Prepackaged CC nutrients in powder form such as that supplied by the
Grand Island Biological Company, 3175
Staley Road, Grand Island, N. Y. 14072
may be used. Earle's lactalbumin hydroly-
sate (ELH), Catalog No. M-11, with
Earle's salts, but without sodium bicarb-
onate is added at the rate of 13.6 gm per
L. The basic formula for ELH is given in
Appendix A.

1.3 Serum: 800 ml (8%) bovine serum
and 200 ml (2%) of fetal
calf serum are first
filtered through cotton
and then added.

1. Report of the Committee on the

1911

1. A report to the

(in the case of the

1. Report of the

1. Report of the

1. Report of the

1. Report of the

1. Report of the

1. Report of the

1. Report of the

1. Report of the

1. Report of the

1. Report of the

1. Report of the

1. Report of the

1. Report of the

1. Report of the

1. Report of the

1. Report of the

1.4 Sodium bicarbonate (NaHCO_3): 22 gm are weighed and added with shaking.

1.5 Vitamins: Difco Laboratories (Detroit, Michigan, USA) TC Vitamins, Eagle, Dried, No. 5879-24 are added (1 gm vial). The basic formula for the vitamins is given in Appendix B.

1.6 Antibiotics: Sodium penicillin G, N.F. No. 8249, the Upjohn Company, Kalamazoo, Michigan 49001, USA: Dihydrostreptomycin sulfate B grade, No. 3021, Calbiochem, Los Angeles, Calif., USA and Mycostatin Sterile Powder (Nystatin, E. R. Squibb and Sons, New York, USA) are added as follows:

Penicillin: To a vial containing a million units of lyophilized sodium penicillin are added 10 ml of $\text{D H}_2\text{O}$; 8 ml of this per batch are added to give a final concentration of 100 units per ml.

Streptomycin: To a 5 gm vial of dihydrostreptomycin sulfate are added 20 ml of D H₂O; 4 ml of this per batch are added to give a final concentration of 100 micrograms per ml.

Mycostatin: To a vial containing 500,000 units, 10 ml of D H₂O are added; 1 ml is used per batch to give a final concentration of 50 units per ml.

1.7 Maintenance medium for seed virus is prepared with 2% fetal bovine serum. (The basic medium, minus serum, may be used as a diluent for the purification of virus.)

1.8 This work should be carried out in hoods, cubicles or other areas free of dust and aerosol contaminants.

2. Filtration: The medium is filtered through a double pass Hormann* filter containing two clarifying pads, followed by filtration through a 4.5 micron and a 2.2 micron Millipore filter in series under 9 lbs. nitrogen pressure. The fluid passes into 10 liter presterilized aspirator bottles (each containing a magnetic stirring bar) with glass filling bells attached.

2.1 All components are presterilized by autoclaving; hoses and connections are flamed with a portable Bunsen burner at time of assembly.

*F. R. Hormann Co., Inc., P. O. Box 229,
Milldale, Conn. 06467.

T. H. M. ...
Box 500

2.2 Frame B-6 shows drums of Baxter bottles on roller mills in a walk-in 37°C incubator. The Vero line of monkey kidney cells is grown in confluent monolayers in these bottles and used to propagate ASF virus for use as antigen in the IEOP and AGDP tests, as well as for other research purposes.

3. Preparation of Vero Cell Cultures

3.1 A saline, trypsin, versene (STV) solution is used for removal of cells from bottles. The formula for the solution is given in Appendix C.


3.2 Trypsinization of Seed Cells:

Thirty ml of STV (in 1 L aspirator bottle with filling bell attached) prewarmed to 37°C is added to each Vero cell culture in a Baxter bottle. When the cells have detached, the suspension is poured into a 2 L aspirator bottle containing 1 L of growth medium as previously described.

1. The first part of the report is devoted to a general survey of the situation in the country. It is a very interesting and well-written account of the country and its people. The author has done a great deal of research and has gathered a wealth of material. The report is a valuable contribution to the knowledge of the country and its people.

2. The second part of the report is devoted to a detailed study of the economic situation. It is a very thorough and well-written account of the economic situation. The author has done a great deal of research and has gathered a wealth of material. The report is a valuable contribution to the knowledge of the economic situation.

3. The third part of the report is devoted to a detailed study of the social situation. It is a very thorough and well-written account of the social situation. The author has done a great deal of research and has gathered a wealth of material. The report is a valuable contribution to the knowledge of the social situation.



3.2 continued

Details of the method are given:

The Baxter bottles are removed from the incubator and taken to a cubicle. All manipulations are carried out with aseptic precautions, using a Bunsen burner to flame stoppers, mouths of containers, etc. Spent medium is poured off and the STV solution is dispensed with a filling bell. After this operation, either the entire drum or separate bottles are rotated to permit the STV solution to reach all cells equally.

From time to time the cells are inspected visually to observe the results of trypsinization and the bottles shaken vigorously to remove the cells.

Partially covered by the weather and given

the weather was not too bad and removed from the

in the afternoon and the weather was not too bad

the weather was not too bad and removed from the

the weather was not too bad and removed from the

the weather was not too bad and removed from the

the weather was not too bad and removed from the

the weather was not too bad and removed from the

the weather was not too bad and removed from the

the weather was not too bad and removed from the

the weather was not too bad and removed from the

in the afternoon

the weather was not too bad and removed from the

the weather was not too bad and removed from the

the weather was not too bad and removed from the

the weather was not too bad and removed from the

3.3 Terminating the Trypsinization:

Chilled ELH medium with 8% bovine and 2% fetal calf serum is placed on a magnetic stirrer. The fluid and trypsinized cells are then added to this from each Baxter bottle by decanting. Presence of serum will terminate the action of the STV.

3.4 Preparation of the New Culture:

Five hundred ml of the decanted mixture containing the trypsinized cells (3.3) is added to 7½ L of growth fluid previously filtered into the 10 L aspirator bottles. The solution is mixed well on a magnetic stirrer and dispensed in 200 ml amounts into sterile empty Baxter bottles held in roller drums. The drums are loaded onto the roller mill at 37°C and allowed to rotate at 10 revolutions per hour for 4 days to achieve confluency.

3. The Commission will report on the progress of the work.

4. Infection of Cell Cultures

4.1 Drums of Baxter bottles with confluent cell cultures are removed from the roller mill to a cubicle. Growth medium is decanted and 10 ml of seed virus at high titer (about $10^{7.8}$ HA_u per ml) added by means of an automatic syringe and 1 L aspirator bottle on a magnetic stirrer.

The drums of Baxter bottles with the infected cells are returned to the incubator for a 4-hour adsorption period. The roller mill is operated at 60 rotations per hour.

4.2 Maintenance medium added. After adsorption of the virus, each bottle receives 10 ml of fresh medium with serum.

4.3 Variations in procedure when CC are used for other types of virus production; (these procedures are not used for antigen production)

-- for stock virus production
each bottle receives 50 ml of medium containing 2% fetal calf serum.

-- for virus purification, the cells are washed with 30 ml of phosphate buffered solution (PBS¹) by gentle rotation. The PBS¹ is poured off and 10 ml of maintenance fluid are added. The formula for PBS¹ is given in Appendix D.

4.3. Weibull in Probability

on test for the time to failure of a component. The Weibull distribution is a generalization of the exponential distribution and is used to model the time to failure of a component.

The Weibull distribution is a two-parameter distribution. The parameters are the shape parameter k and the scale parameter λ . The probability density function (PDF) of the Weibull distribution is given by:

The cumulative distribution function (CDF) of the Weibull distribution is given by:

The Weibull distribution is used to model the time to failure of a component. The shape parameter k determines the shape of the distribution. The scale parameter λ determines the scale of the distribution. The Weibull distribution is a generalization of the exponential distribution.

5. Harvesting the ASF Virus

5.1 The harvesting method is simple. For the usual IEOP and AGDP antigen or for viral seed material, the cells are harvested with the fluid when the cells show cytopathic effect (CPE) and are sloughing from the walls of the bottles. This usually occurs between 48 and 72 hours after inoculation.

5.2 Where purer viral antigen or virus is required, the fluid portion only in the Baxter bottles is harvested at 24 hours postinoculation.

2. How is the ...

2.1. The ... is ...

For the ...

but also ...

has ...

the ...

along ...

it ...

to ...

2.2. When ...

when ...

only ...

at ...

5.3 Harvesting Technique: Each bottle is shaken vigorously to remove remaining cells and the cell-fluid mixture decanted. The mixture is then centrifuged at 450 X G in the cold (6°C) for 30 minutes. A centrifuge* with an 8 inch radius measured to 1/2 the depth of fluid in the centrifuge tube was used at a speed of 2,000 revolutions per minute (r.p.m.) at 6°C. (Appendix E has a graph for converting relative centrifugal force $\overline{R.c.f.}$ or \overline{G} to r.p.m. or vice versa.)

The pellet, which contains the cellular debris and most of the virus, is saved for further treatment and use as IEOP and AGDP antigen. (The fluid portion also contains some virus and may be retained, but it is not used for antigen.)

*Model PRJ of the International Equipment Company, Needham, Mass., USA.













6. Completion of Antigen Preparation

6.1 The packed cell debris is resuspended in 2 volumes of PBS and sonicated for 2 minutes (C1,C2). This is followed by ultracentrifugation in a Spinco* 40 rotor for 60 minutes at 35,000 r.p.m. (C3,C4). The supernate is carefully removed from the ultracentrifuge tube by syringe and canula (C5,C6) and transferred to another container for storage at -70°C . Sodium azide (Na_3N) is added to give a final concentration of 0.1%.

* Beckman Instruments, Inc., Spinco Division, Route 22 at Summit Boulevard, Mountainside, New Jersey 07091.

The first of these is the fact that the
 results of the experiments are in good
 agreement with the theoretical predictions.
 The second is the fact that the
 results are in good agreement with the
 results of the experiments of other
 workers. The third is the fact that
 the results are in good agreement with
 the results of the experiments of other
 workers. The fourth is the fact that
 the results are in good agreement with
 the results of the experiments of other
 workers. The fifth is the fact that
 the results are in good agreement with
 the results of the experiments of other
 workers. The sixth is the fact that
 the results are in good agreement with
 the results of the experiments of other
 workers. The seventh is the fact that
 the results are in good agreement with
 the results of the experiments of other
 workers. The eighth is the fact that
 the results are in good agreement with
 the results of the experiments of other
 workers. The ninth is the fact that
 the results are in good agreement with
 the results of the experiments of other
 workers. The tenth is the fact that
 the results are in good agreement with
 the results of the experiments of other
 workers.

The first of these is the fact that the
 results of the experiments are in good
 agreement with the theoretical predictions.
 The second is the fact that the
 results are in good agreement with the
 results of the experiments of other
 workers. The third is the fact that
 the results are in good agreement with
 the results of the experiments of other
 workers. The fourth is the fact that
 the results are in good agreement with
 the results of the experiments of other
 workers. The fifth is the fact that
 the results are in good agreement with
 the results of the experiments of other
 workers. The sixth is the fact that
 the results are in good agreement with
 the results of the experiments of other
 workers. The seventh is the fact that
 the results are in good agreement with
 the results of the experiments of other
 workers. The eighth is the fact that
 the results are in good agreement with
 the results of the experiments of other
 workers. The ninth is the fact that
 the results are in good agreement with
 the results of the experiments of other
 workers. The tenth is the fact that
 the results are in good agreement with
 the results of the experiments of other
 workers.

THIS COMPLETES THE FIRST STEP -- THAT OF
ANTIGEN PREPARATION -- FOR EITHER THE
IEOP OR AGDP TESTS.

* * * * *

NEXT, PROTOCOLS FOR PERFORMING THE IEOP
TESTS WILL BE GIVEN.

THE UNIVERSITY OF CHICAGO
LIBRARY
540 EAST 57TH STREET
CHICAGO, ILL. 60637

* 4 - 1 - 1

THE UNIVERSITY OF CHICAGO
LIBRARY
540 EAST 57TH STREET
CHICAGO, ILL. 60637

7. Performing the IEOP Test

7.1 Equipment

Any suitable immunoelectrophoresis apparatus with a power pack developing 500 volts may be used. The tests described here were carried out with a Gelman Power Supply No. 38206 (B-8) and accessories for immunoelectrophoresis (Basic Outfit No. 51841 and Screening Kit No. 51482 are typical).

7.2 Buffer Solution

An electrophoresis (barbital) buffer solution of 0.1 ionic strength, pH 8.6 is prepared according to the formula given in Appendix F.

7.3 Agar Gel

An agarose or Noble (Difco, Detroit, Michigan, USA) agar is prepared for coating the slides according to the formula given in Appendix F.

V. THE CONSTITUTION

1.1. Introduction

The Constitution of the United States is the supreme law of the land. It is the foundation of the government and the rights of the people. The Constitution is a living document that has been interpreted and amended over time. It is the source of the powers of the federal government and the states. It is the source of the rights of the individual. It is the source of the structure of the government. It is the source of the relationship between the federal government and the states. It is the source of the relationship between the government and the people.

1.2. The Preamble

We the People of the United States, in Order to form a more perfect Union, establish Justice, insure domestic Tranquility, provide for the common defence, promote the general Welfare, and secure the Blessings of Liberty to ourselves and our Posterity, do hereby adopt this Constitution.

1.3. The Principles of the Constitution

The Constitution is based on the principles of federalism, separation of powers, and checks and balances. Federalism is the division of power between the federal government and the states. Separation of powers is the division of power between the executive, legislative, and judicial branches. Checks and balances are the mechanisms by which each branch of government can limit the power of the other branches.

1.4. The Bill of Rights

The Bill of Rights is the first ten amendments to the Constitution. It was added to the Constitution in 1791. It guarantees the rights of the individual against the federal government.

1.5. The Fourteenth Amendment

The Fourteenth Amendment is the first amendment to the Bill of Rights. It was added to the Constitution in 1868. It guarantees the rights of the individual against the states.

The Fourteenth Amendment is the first amendment to the Bill of Rights. It was added to the Constitution in 1868. It guarantees the rights of the individual against the states.

The Fourteenth Amendment is the first amendment to the Bill of Rights. It was added to the Constitution in 1868. It guarantees the rights of the individual against the states.

7.4 COMMENT ON FRAMES C-7, C-8 and C-9

FRAME C-11 shows the chief mechanical devices used in the IEOP test. To the left rear is the POWER CONVERTER or POWER SUPPLY. In front of it is a LEVELING TABLE SET with a SPIRIT LEVEL in the center of the table. On the right front is an ELECTROPHORESIS CHAMBER.

FRAME C-12 shows another model of the power converter and electrophoresis chamber.

FRAME D-1 is a highly diagrammatic sketch of the power converter and electrophoresis chamber to aid in identification of the components and use of the units.

THEY ARE ALWAYS NO MORE

1

THEY ARE ALWAYS NO MORE

THEY ARE ALWAYS NO MORE

THEY ARE ALWAYS NO MORE

THEY ARE ALWAYS NO MORE

THEY ARE ALWAYS NO MORE

THEY ARE ALWAYS NO MORE

THEY ARE ALWAYS NO MORE

THEY ARE ALWAYS NO MORE

THEY ARE ALWAYS NO MORE

THEY ARE ALWAYS NO MORE

THEY ARE ALWAYS NO MORE


THEY ARE ALWAYS NO MORE


THEY ARE ALWAYS NO MORE

THEY ARE ALWAYS NO MORE


THEY







POWER CONVERTER




SLIDE RACK

1932



1932



1932


1932

1932

8. Preparing the Slides


With the Gelman apparatus shown, in Frames C-11 and C-12, a maximum of 192 serum samples may be tested simultaneously. Six ordinary glass slides (about 2.5 cm by 7.5 cm) each are held in a frame divided into two sections with three slides placed end to end in each section. Two movable front screws of the leveling table are adjusted until the bubble of the removable spirit level is centered. The frame is then firmly clamped on the leveling table. This assures even thickness of agar gel on the glass slides.

the first of these is the fact that the
the second is the fact that the
the third is the fact that the
the fourth is the fact that the
the fifth is the fact that the
the sixth is the fact that the
the seventh is the fact that the
the eighth is the fact that the
the ninth is the fact that the
the tenth is the fact that the
the eleventh is the fact that the
the twelfth is the fact that the
the thirteenth is the fact that the
the fourteenth is the fact that the
the fifteenth is the fact that the
the sixteenth is the fact that the
the seventeenth is the fact that the
the eighteenth is the fact that the
the nineteenth is the fact that the
the twentieth is the fact that the
the twenty-first is the fact that the
the twenty-second is the fact that the
the twenty-third is the fact that the
the twenty-fourth is the fact that the
the twenty-fifth is the fact that the
the twenty-sixth is the fact that the
the twenty-seventh is the fact that the
the twenty-eighth is the fact that the
the twenty-ninth is the fact that the
the thirtieth is the fact that the
the thirty-first is the fact that the
the thirty-second is the fact that the
the thirty-third is the fact that the
the thirty-fourth is the fact that the
the thirty-fifth is the fact that the
the thirty-sixth is the fact that the
the thirty-seventh is the fact that the
the thirty-eighth is the fact that the
the thirty-ninth is the fact that the
the fortieth is the fact that the
the forty-first is the fact that the
the forty-second is the fact that the
the forty-third is the fact that the
the forty-fourth is the fact that the
the forty-fifth is the fact that the
the forty-sixth is the fact that the
the forty-seventh is the fact that the
the forty-eighth is the fact that the
the forty-ninth is the fact that the
the fiftieth is the fact that the
the fifty-first is the fact that the
the fifty-second is the fact that the
the fifty-third is the fact that the
the fifty-fourth is the fact that the
the fifty-fifth is the fact that the
the fifty-sixth is the fact that the
the fifty-seventh is the fact that the
the fifty-eighth is the fact that the
the fifty-ninth is the fact that the
the sixtieth is the fact that the
the sixty-first is the fact that the
the sixty-second is the fact that the
the sixty-third is the fact that the
the sixty-fourth is the fact that the
the sixty-fifth is the fact that the
the sixty-sixth is the fact that the
the sixty-seventh is the fact that the
the sixty-eighth is the fact that the
the sixty-ninth is the fact that the
the seventieth is the fact that the
the seventy-first is the fact that the
the seventy-second is the fact that the
the seventy-third is the fact that the
the seventy-fourth is the fact that the
the seventy-fifth is the fact that the
the seventy-sixth is the fact that the
the seventy-seventh is the fact that the
the seventy-eighth is the fact that the
the seventy-ninth is the fact that the
the eightieth is the fact that the
the eighty-first is the fact that the
the eighty-second is the fact that the
the eighty-third is the fact that the
the eighty-fourth is the fact that the
the eighty-fifth is the fact that the
the eighty-sixth is the fact that the
the eighty-seventh is the fact that the
the eighty-eighth is the fact that the
the eighty-ninth is the fact that the
the ninetieth is the fact that the
the ninety-first is the fact that the
the ninety-second is the fact that the
the ninety-third is the fact that the
the ninety-fourth is the fact that the
the ninety-fifth is the fact that the
the ninety-sixth is the fact that the
the ninety-seventh is the fact that the
the ninety-eighth is the fact that the
the ninety-ninth is the fact that the
the hundredth is the fact that the



A small amount of melted agar gel is introduced at the junctures of the slides and allowed to spread by capillary action between the slides and the frame. The agar gel is allowed to solidify and form a seal which will prevent leakage when the agar gel is applied to the surface. Ten ml of melted agar gel are then pipetted on each section of three slides to form a continuous layer of uniform thickness covering the three slides. This is repeated on the other three slides of the frame. The agar gel is allowed to solidify for 30 minutes before use.





9. Preparing the Wells

Although wells may be cut by cork borers, using a pattern, greater accuracy is obtained when a commercially available pattern cutter is employed. The wells are 3 mm in diameter and spaced at 10 mm between edges of wells. The gel plugs are removed by a 12 gauge canula, or Pasteur pipette connected with tubing to a vacuum source. A trap to prevent the agar gel plus from entering the vacuum source hose is made by interposing a flask between the canula hose and the vacuum source line.

the water level in the
lake was about 10 feet
higher than in the
river. The water in the
lake was very muddy and
the water in the river
was very clear. The water
in the lake was very
warm and the water in
the river was very cold.
The water in the lake
was very muddy and the
water in the river was
very clear. The water
in the lake was very
warm and the water in
the river was very cold.
The water in the lake
was very muddy and the
water in the river was
very clear. The water
in the lake was very
warm and the water in
the river was very cold.






10. Preparing Special Filling Pipettes

A special pipette is used to insure greater speed and accuracy in filling the wells. Ordinary Pasteur pipettes are rapidly converted to angle-tip micropipettes. The tip of the regular pipette is heated above a Bunsen burner, pulled out to about 6 inches and the glass bent to an approximate 45° angle. It is raised to cool the glass and the tip bent inward until the microtube portion breaks; usually the bent portion remaining is about 2 cm in length.








11. Filling the Wells

The wells may be filled accurately and rapidly with the special micro-pipette and a bulb (E-1). Note that for extra steadiness, the hand may be rested on the bench surface. A separate pipette must be used for each sample. (The antigen wells may all be filled with one pipette.)

Suspect ASF antiserum is placed in one of a pair of wells nearest the anode (+). The antigen is placed in the other well of the pair nearest the cathode (-).



11.1 COMMENT ON FRAME E-4

The Gelman power supply and electrophoresis chamber are shown in the photograph. One immuno-frame is in place within the chamber, which will hold four. The manufacturer recommends that not more than three frames be used at a time in some chambers, such as the one shown in Frame C-12. Comparison with the drawing in the next frame after the photograph will make components and their use clear. (Please note that the drawing is somewhat simplified and diagrammatic; for instance, the wells on each side are more extensive than shown.) The placement of the double wicks is best shown in the drawing. Only the top of the wicks, white bands across each end of the immuno-frame, can be seen in the photograph.

6. 1994 1995 1996 1997 1998 1999 2000 2001 2002 2003 2004 2005 2006 2007 2008 2009 2010 2011 2012 2013 2014 2015 2016 2017 2018 2019 2020 2021 2022 2023 2024 2025 2026 2027 2028 2029 2030 2031 2032 2033 2034 2035 2036 2037 2038 2039 2040 2041 2042 2043 2044 2045 2046 2047 2048 2049 2050 2051 2052 2053 2054 2055 2056 2057 2058 2059 2060 2061 2062 2063 2064 2065 2066 2067 2068 2069 2070 2071 2072 2073 2074 2075 2076 2077 2078 2079 2080 2081 2082 2083 2084 2085 2086 2087 2088 2089 2090 2091 2092 2093 2094 2095 2096 2097 2098 2099 2100 2101 2102 2103 2104 2105 2106 2107 2108 2109 2110 2111 2112 2113 2114 2115 2116 2117 2118 2119 2120 2121 2122 2123 2124 2125 2126 2127 2128 2129 2130 2131 2132 2133 2134 2135 2136 2137 2138 2139 2140 2141 2142 2143 2144 2145 2146 2147 2148 2149 2150 2151 2152 2153 2154 2155 2156 2157 2158 2159 2160 2161 2162 2163 2164 2165 2166 2167 2168 2169 2170 2171 2172 2173 2174 2175 2176 2177 2178 2179 2180 2181 2182 2183 2184 2185 2186 2187 2188 2189 2190 2191 2192 2193 2194 2195 2196 2197 2198 2199 2200 2201 2202 2203 2204 2205 2206 2207 2208 2209 2210 2211 2212 2213 2214 2215 2216 2217 2218 2219 2220 2221 2222 2223 2224 2225 2226 2227 2228 2229 2230 2231 2232 2233 2234 2235 2236 2237 2238 2239 2240 2241 2242 2243 2244 2245 2246 2247 2248 2249 2250 2251 2252 2253 2254 2255 2256 2257 2258 2259 2260 2261 2262 2263 2264 2265 2266 2267 2268 2269 2270 2271 2272 2273 2274 2275 2276 2277 2278 2279 2280 2281 2282 2283 2284 2285 2286 2287 2288 2289 2290 2291 2292 2293 2294 2295 2296 2297 2298 2299 2300 2301 2302 2303 2304 2305 2306 2307 2308 2309 2310 2311 2312 2313 2314 2315 2316 2317 2318 2319 2320 2321 2322 2323 2324 2325 2326 2327 2328 2329 2330 2331 2332 2333 2334 2335 2336 2337 2338 2339 2340 2341 2342 2343 2344 2345 2346 2347 2348 2349 2350 2351 2352 2353 2354 2355 2356 2357 2358 2359 2360 2361 2362 2363 2364 2365 2366 2367 2368 2369 2370 2371 2372 2373 2374 2375 2376 2377 2378 2379 2380 2381 2382 2383 2384 2385 2386 2387 2388 2389 2390 2391 2392 2393 2394 2395 2396 2397 2398 2399 2400 2401 2402 2403 2404 2405 2406 2407 2408 2409 2410 2411 2412 2413 2414 2415 2416 2417 2418 2419 2420 2421 2422 2423 2424 2425 2426 2427 2428 2429 2430 2431 2432 2433 2434 2435 2436 2437 2438 2439 2440 2441 2442 2443 2444 2445 2446 2447 2448 2449 2450 2451 2452 2453 2454 2455 2456 2457 2458 2459 2460 2461 2462 2463 2464 2465 2466 2467 2468 2469 2470 2471 2472 2473 2474 2475 2476 2477 2478 2479 2480 2481 2482 2483 2484 2485 2486 2487 2488 2489 2490 2491 2492 2493 2494 2495 2496 2497 2498 2499 2500 2501 2502 2503 2504 2505 2506 2507 2508 2509 2510 2511 2512 2513 2514 2515 2516 2517 2518 2519 2520 2521 2522 2523 2524 2525 2526 2527 2528 2529 2530 2531 2532 2533 2534 2535 2536 2537 2538 2539 2540 2541 2542 2543 2544 2545 2546 2547 2548 2549 2550 2551 2552 2553 2554 2555 2556 2557 2558 2559 2560 2561 2562 2563 2564 2565 2566 2567 2568 2569 2570 2571 2572 2573 2574 2575 2576 2577 2578 2579 2580 2581 2582 2583 2584 2585 2586 2587 2588 2589 2590 2591 2592 2593 2594 2595 2596 2597 2598 2599 2600 2601 2602 2603 2604 2605 2606 2607 2608 2609 2610 2611 2612 2613 2614 2615 2616 2617 2618 2619 2620 2621 2622 2623 2624 2625 2626 2627 2628 2629 2630 2631 2632 2633 2634 2635 2636 2637 2638 2639 2640 2641 2642 2643 2644 2645 2646 2647 2648 2649 2650 2651 2652 2653 2654 2655 2656 2657 2658 2659 2660 2661 2662 2663 2664 2665 2666 2667 2668 2669 2670 2671 2672 2673 2674 2675 2676 2677 2678 2679 2680 2681 2682 2683 2684 2685 2686 2687 2688 2689 2690 2691 2692 2693 2694 2695 2696 2697 2698 2699 2700 2701 2702 2703 2704 2705 2706 2707 2708 2709 2710 2711 2712 2713 2714 2715 2716 2717 2718 2719 2720 2721 2722 2723 2724 2725 2726 2727 2728 2729 2730 2731 2732 2733 2734 2735 2736 2737 2738 2739 2740 2741 2742 2743 2744 2745 2746 2747 2748 2749 2750 2751 2752 2753 2754 2755 2756 2757 2758 2759 2760 2761 2762 2763 2764 2765 2766 2767 2768 2769 2770 2771 2772 2773 2774 2775 2776 2777 2778 2779 2780 2781 2782 2783 2784 2785 2786 2787 2788 2789 2790 2791 2792 2793 2794 2795 2796 2797 2798 2799 2800 2801 2802 2803 2804 2805 2806 2807 2808 2809 2810 2811 2

and it would not be possible to do so.

[illegible]

act before the day when the ...

9109 Jan 1971 Brown 2000 1000 1000 1000

Journal of Management Studies, 19(6), 701-718.



12. Electrophoresis

The bridges of the electrophoresis chamber should be adjusted to the second position from the outer edges of the chamber. Four frames of 6 slides each may be placed across the support bridges, although only one is shown in the illustration. Wicks (for conduction of current through the gel) of cellulose acetate are saturated in the electrophoresis buffer solution described. One end of a wick is placed on each end of the agar and the other end of the wick allowed to rest in the buffer solution of the chamber. Double wicking or placing two thicknesses of wick is recommended.

10. The first of these is the

fact that the effect of the

operation is to adjust the

position of the other end of

the bar. The amount of this

adjustment is the amount of the

change in the position of the

other end of the bar.

The amount of this change

is the amount of the

change in the position of the

other end of the bar.

The amount of this change

is the amount of the

change in the position of the

other end of the bar.

11.

The usual voltage setting for electrophoresis is 19 volts per cm or 450 volts when using the Gelman apparatus. With some chambers this voltage must be reduced when three immuno-frames are used. The test can usually be performed satisfactorily with voltages between 300 and 450. The most satisfactory length of time at PIADL has been 30 minutes. To determine the optimal duration for particular apparatus and reagents, test runs of from 15 to 60 minutes may be made.

The life of a buffer solution may be extended by alternating the anode and cathode. With the Gelman apparatus this is done simply by the setting of a knob. With other types of apparatus lacking such a device, the chamber might be reversed. When the direction of electric current is thus changed, the wells must be filled in the opposite manner, making sure that the serum is on the side nearest the anode and the antigen on the side nearest the cathode. By reversing the current flow in this manner the buffer solution may be used for 10 or more electrophoretic runs. When a total of 4 frames are used each time, the longevity of the buffer is increased proportionately. For example, if only 1 frame were used for each run, the life of the buffer would be increased fourfold.

13. Final Checks Before Electrophoresis

The antisera must be in the wells nearest the anode (+) and the antigen nearest the cathode (-). The direction of electric current flow is determined to be in the appropriate direction. No chamber components should be touched while the current is flowing as there is danger of electric shock. (The Gelman apparatus is equipped with a safety feature which turns off the electric current to the chamber when the lid is lifted.)

At the end of electrophoresis, the current is shut off and the frames removed for viewing by indirect transmitted light. A viewing box is convenient; the one shown has an adjustable mirror. A precipitin line or lines between the two wells indicates a positive test.

1. The first of these is the

fact that the world is

not a uniform whole

but a collection of

many different parts

which are not

connected by any

common principle

of unity or order

and which are

not subject to

any common law

or common rule

of conduct

or common principle

of action or

of thought or

of feeling or

of any other

kind of

activity or

of any other

SORVAL
SUPERSPEED
RC2-B




Lisbon

Ag



Rome

Ag



+

-

For a permanent record, the wet agar gel slide may be photographed. Slides may also be preserved by washing for a 48-hour-period each time, in two changes of 2% NaCl solution, followed by a 1-hour rinse in D H₂O. Slides are then air dried and stained with amido (Amido Black 10B dye, Bayerwerke, Leverkusen, Germany) or buffalo (Buffalo Black, NBR, Allied Chemicals, New York) black. Six gm of one of these dyes are dissolved in a mixture of methyl alcohol (450 ml), D H₂O (450 ml) and glacial acetic acid (100 ml); then filtered through two layers of gauze before use. The decolorizing solution is made by the same formula with omission of the dye. Decolorizing is carried out until the agar gel background is clear and the lines well delineated. Slides are subsequently numbered directly on the

dried agar film with India ink and stored in a dry place. This slide can be used as a "negative" in a photographic enlarger and the image printed on photographic paper. The negative image resembles the lines as viewed by indirect light before staining; these reproductions are suitable for record or for publication.

... and ...
...
...
...
...
...
...
...
...
...

This completes the description of the IEOP test.

* * * * *

The next section of this microfiche gives protocols for performing the AGDP test. This test is similar to the IEOP test, but is much simpler and requires little equipment and no special laboratory items. The same antigen is used for both tests. The next frame summarizes the basic steps in the AGDP test.

of the month of the year 1897

1897

1897

The first of the month of the year 1897
the first of the month of the year 1897
the first of the month of the year 1897
the first of the month of the year 1897
the first of the month of the year 1897
the first of the month of the year 1897
the first of the month of the year 1897
the first of the month of the year 1897

The Agar Gel Diffusion Precipitation Test

The AGDP test is useful as a research and diagnostic tool. If adequate amounts of either antigen or antisera are available in the fluids or tissues of the suspect animal, the test may be performed in from 18 to 48 hours. The technique is simple and requires no elaborate equipment or costly supplies. If a known ASF antisera is available, the test may be performed by placing fluids or tissues of an animal which died of suspected ASF in wells opposite the antisera. Viral antigen has been demonstrated from the blood, spleen, liver, and kidney of such animals.

If a known ASF antigen, such as that made for the IEOP test, is available, sera or other fluids from suspect animals may be used to detect specific antibodies.

The Age of the Earth and the History of the World

CHAPTER I. OF THE ORIGIN OF THE EARTH

THE ORIGIN OF THE EARTH is a subject which has long attracted the attention of philosophers and naturalists. The various hypotheses which have been proposed to explain the origin of the world, and the manner in which it has been formed, are all founded on conjecture and speculation. The most plausible of these hypotheses is that which supposes the earth to have been formed out of a mass of matter, which was originally in a state of confusion and disorder. This mass of matter, it is supposed, was first set in motion by some powerful agent, and then gradually condensed into the form of the earth.

THE ORIGIN OF THE EARTH is a subject which has long attracted the attention of philosophers and naturalists. The various hypotheses which have been proposed to explain the origin of the world, and the manner in which it has been formed, are all founded on conjecture and speculation. The most plausible of these hypotheses is that which supposes the earth to have been formed out of a mass of matter, which was originally in a state of confusion and disorder. This mass of matter, it is supposed, was first set in motion by some powerful agent, and then gradually condensed into the form of the earth.

THE ORIGIN OF THE EARTH is a subject which has long attracted the attention of philosophers and naturalists. The various hypotheses which have been proposed to explain the origin of the world, and the manner in which it has been formed, are all founded on conjecture and speculation. The most plausible of these hypotheses is that which supposes the earth to have been formed out of a mass of matter, which was originally in a state of confusion and disorder. This mass of matter, it is supposed, was first set in motion by some powerful agent, and then gradually condensed into the form of the earth.

The same method of washing slides in saline and D H₂O given in the previous section must always be used, as putrid or macerated tissue may result in non-specific lines when AGDP plates are viewed before adequate washing.

The same result was obtained

using and H₂O₂ added to the solution

and the most common method is to

or reacted with the solution in a

specific type of apparatus

which before use is washed

The agar gel and buffer solution used are identical to those described in detail for the IEOP test. The method of preparing slides is also the same, except that individual slides or Petri dishes are poured without using the immunoframes. Three ml of melted agar per slide or 15 ml of melted agar per 8.5 cm Petri dish are used.

Wells may be cut in a variety of patterns so that antigen and antisera are placed in opposite wells. Frequently a circular pattern is used with the chief reagent placed in the central well and test sera, tissue or fluids in the peripheral wells. Thus the central well is filled with antigen when testing for the presence of antibodies in the peripheral wells. Frame F-7 is a photograph of an AGDP test; in this case a template was

The first part of the report is devoted to a description of the experimental apparatus and the method of observation. The second part contains the results of the experiments and a discussion of the factors which influence the rate of reaction. The third part is devoted to a comparison of the results obtained with those obtained by other investigators. The fourth part contains the conclusions of the author.

The experimental apparatus consists of a glass vessel of known volume, in which a known amount of gas is introduced. The gas is then ignited by a small flame, and the rate of reaction is observed by measuring the time which elapses between the ignition and the extinction of the flame. The results of the experiments show that the rate of reaction is influenced by the pressure, the temperature, and the nature of the gas.

The rate of reaction increases with increasing pressure and temperature. It also increases with increasing molecular weight of the gas. The results obtained by other investigators are in good agreement with those obtained in this experiment.

The conclusions of the author are that the rate of reaction is influenced by the pressure, the temperature, and the nature of the gas. The rate of reaction increases with increasing pressure and temperature, and with increasing molecular weight of the gas.

used to cut the central well and 6 peripheral wells at one time.

The well cutting template consisted of a metal disc to which six 4 mm sharpened metal tubes were affixed in a circle around a central tube of the same diameter. The approximate distance from the outer edge of the central well to the near edge of each peripheral well was 5 mm.

COMMENT ON THE NEXT FRAME:

This is a photograph of a typical AGDP test. Wells were filled as follows:

Center Well: ASFV Antigen

Peripheral Wells 1, 3 and 5: Normal
Swine Serum

Peripheral Wells 2, 4 and 6: Sera
from swine which recovered from
ASF.

The first of these is the fact that the
 number of cases of the disease has
 been increasing steadily since 1910.
 The second is the fact that the disease
 is now found in all parts of the world.
 The third is the fact that the disease
 is now found in all parts of the world.
 The fourth is the fact that the disease
 is now found in all parts of the world.
 The fifth is the fact that the disease
 is now found in all parts of the world.
 The sixth is the fact that the disease
 is now found in all parts of the world.
 The seventh is the fact that the disease
 is now found in all parts of the world.
 The eighth is the fact that the disease
 is now found in all parts of the world.
 The ninth is the fact that the disease
 is now found in all parts of the world.
 The tenth is the fact that the disease
 is now found in all parts of the world.



As in the IEOP test, the agar must be washed properly in saline to remove nonspecific lines. By placing antisera in the central well it is possible to test for the presence of antigen by placing bits of tissue or body fluids from suspect animals in the outer wells. Dr. W. R. Hess of PIADL states that in case of tissues containing ASF virus, it is probable that only one line may appear. For a more complete discussion of the IEOP and AGDP tests for ASF, reference should be made to Pan, De Boer and Hess, 1972 (full reference given in Frames A-1 and A-2).

This completes the description of the
AGDP test.

* * * * *

The next section (Fiche 2) comprises a
brief description of the FA test. This
technique is quite different from the
two previously described. Its basic
purpose is to demonstrate ASF virus or
viral antigen in tissues. The first
frame in the second fiche summarizes
the basic steps in performing this test.

THE UNIVERSITY OF CHICAGO PRESS

CHICAGO, ILL.

1910

THE UNIVERSITY OF CHICAGO PRESS
PUBLISHED BY THE UNIVERSITY OF CHICAGO PRESS
CHICAGO, ILL.
1910
THE UNIVERSITY OF CHICAGO PRESS
PUBLISHED BY THE UNIVERSITY OF CHICAGO PRESS
CHICAGO, ILL.
1910

Appendix A

The formula for Earle's salts in the ELH medium is:

Component	mg/L
NaCl	8,000.0
KCl	400.0
NaH ₂ PO ₄	60.0
CaCl ₂ (anhydrous)	140.0
MgCl ₂ · 6 H ₂ O	100.0
MgSO ₄ · 7 H ₂ O	100.0
KH ₂ PO ₄ (anhydrous)	60.0
Glucose	1,000.0
Lactalbumin Hydrolysate	6,500.0
Phenol red	10.0
NaHCO ₃	350.0

Appendix B

The formula for the amount of added vitamins is:

<u>Vitamin</u>	<u>Amount</u>
D Biotin	0.24 mg
Folic acid	0.44 mg
Niacinamide	0.12 mg
Calcium pantothenic	0.24 mg
Pyridoxal hydrochloride	0.20 mg
Thiaminde hydrochloride	0.34 mg
Riboflavin	0.04 mg
Choline chloride	0.14 mg

Appendix C

Formula for Saline, Trypsin, Versene (STV) Solution

Component	Amount
NaCl	8.0 gm
KCl	0.4 gm
Dextrose	1.0 gm
NaHCO ₃	0.58 gm
Trypsin (Difco 1:250)	0.5 gm
Versene (EPTA) (Disodium salt 1:250)	0.2 gm
D H ₂ O to make	1,000 ml

Filter and store below -20°C; do not
refreeze and use again.

Appendix D


Formula,

Phosphate buffered solution (PBS¹)
for purification of ASF virus

NaCl	8.0	gm
KCl	0.2	gm
Na ₂ HPO ₄ · 7 H ₂ O	2.172	gm
KN ₂ PO ₄	0.2	gm
D H ₂ O	1,000	ml

Appendix E

LABORATORY AID FOR OBTAINING RELATIVE CENTRIFUGAL FORCE (G) or REVOLUTIONS PER MINUTE (RPM)



Appendix F

Electrophoresis Buffer Solution:

(Barbital buffer solution, ionic strength 0.1, pH 8.6).

Sodium barbital	13.38 gm
Sodium acetate (3 H ₂ O)	8.83 gm
Distilled H ₂ O to make	1.50 liters
Adjust pH with conc. HCl to pH 8.6	

Electrophoresis Agar Gel:

6 gm agarose, or 10 gm Noble (Difco, Detroit, Mich., USA) special agar
10 ml, 10% sodium azide
250 ml barbital buffer, above
750 ml distilled water
Heat to 100°C until agarose is melted for use.

